

# Pet Cooling Gel Bandana

## Product function / characteristics:

1. MATERIAL : with breathable mesh instant cooling fabric , soft and comfortable . Contains of ice crystal power inside .
2. Adjustable sizes for fitting different dog breeds

## Use tips :

Just soak the bandana in cold water and wring out, and the cooling crystals activate, staying hydrated for several days.

You can also put into refrigerator to get better cooling effect .

The High Tech Instant cooling fabric provides comfort around your dog's neck, while the super-absorbent cooling crystals retain moisture to keep the wearer cool all day long.

# Pet Cooling Gel Neck Band



COLDING



COMFORTABLE



BREATHABLE




DRIPPING WET

Keep the wearer cool all day long

High Tech

Polyester fiber

Contains cooling gel powder



# Instructions



1. Immerse Gel Neck in water.



2. After soaking for 2-3 minutes, remove and gently twist the excess water.



3. It's better to use it after refrigerated for 10 minutes.



4. Keep the wearer cool all day long.

## Where to use?

Indoor



Walk a dog



Playing



Climbing





## Product information

**Brand:** DogLemi

**Name:** Pet Cooling Gel Neck Band

**Item NO:** PD60110

**Size:** S/M/L

**Color:** Blue

**Material:** Polyester fiber

### ▼ SIZE CHART



### SIZE CHART

SIZE	LENGTH	WIDTH
S	25.2" / 64cm	2" / 5cm
M	31.5" / 80cm	2.4" / 6cm
L	38.6" / 98cm	2.8" / 7cm



## Model display







